

Bronchial Stents with Rings

Code No:	A O.D.(mm)	B I.D.(mm)	C Length (mm)	D Flanged O.D. (mm)
BS-0613-S	6	5	13	8
BS-0813-S	8	5.5	13	10
BS-1013-S	10	8	13	13
BS-0619-S	6	5	19	8
BS-0819-S	8	5.5	19	10
BS-1019-S	10	8	19	13
BS-0313-S	Length 13, 1 each of O.D. 6, 8, 10mm			
BS-0319-S	Length 19, 1 each of O.D. 6, 8, 10mm			

