

Radiopaque Standard Tracheal T-Tubes

Code No:	O.D. (mm)	A	B	C	D
RST-06-S	6	10	6	48	5
RST-07-S	7	13	9	48	6
RST-08-S	8	16	12	36	6
RST-09-S	9	19	15	40	8
RST-10-S	10	23	17	44	8
RST-11-S	11	27	20	50	9
RST-12-S	12	31	21	55	11
RST-13-S	13	31	22	60	11
RST-14-S	14	32	23	64	11
RST-15-S	15	34	24	64	11
RST-16-S	16	31	23	63	12

Radiopaque Long Tracheal T-Tubes

Code No:	O.D. (mm)	A	B	C	D
RLT-06-S	6	35	10	38	5
RLT-07-S	7	42	10	38	6
RLT-08-S	8	60	13	31	6.5
RLT-09-S	9	65	15	34	8
RLT-10-S	10	70	18	41	8
RLT-11-S	11	60	30	61	9
RLT-12-S	12	75	22	49	11
RLT-13-S	13	75	23	59	11
RLT-14-S	14	75	25	59	11
RLT-16-S	16	85	25	59	12
RLT-18-S	18	85	25	59	14

Radiopaque Thoracic Tracheal T-Tubes (extra-long)

Code No:	O.D. (mm)	A	B	C	D
RELT-08-S	8	61	50	40	6.5
RELT-10-S	10	70	50	41	8
RELT-12-S	12	100	40	64	11
RELT-14-S	14	99	50	74	11