

Neonatal and Pediatric Tracheostomy Tubes

No one does it better!

A complete line of Tracheostomy tubes for newborns and children is an indispensable part of any practice. Due to the specific and delicate anatomy of these little patients, specific design requirements and materials are required — No one does it better!



Each tube is supplied with a wide neck strap.

- RING-SHAPED handle of the obturator which ensures safe guidance during the insertion of the tube.
- Back of hub cut off at a 45° for better patient comfort.
- Soft, compliant and flexible, “Less Atraumatic.” Reduces the possibility of pressure sores.
- Wings of neck flange form a 150° angle towards the head. Ideally adapted to the specific anatomical requirements of neonates and children, resulting in a more comfortable fit.
- The tubes are individually packaged with a neck strap tailor-made for neonates and children, making changing trach’s less cumbersome and easier for families.
- Thinner walls than most other tubes. Optimal ratio between tube size and lumen. This ratio is most important for very small diameters since the resistance to air flow increases dramatically.
- Radiopaque can remain in situ during X-ray monitoring; X-ray pictures will show the entire tube — not only an X-ray contrast.
- Neck flange is supplied with large eyelets for patient comfort and ease of changing neck ties.

*Used and endorsed by
Cincinnati Children’s Otolaryngology Head and Neck Surgery Dept.*



Bryan Medical, Inc. (513) 272-1600 (bus) (513) 272-1610 (fax)
3914 Miami Road, Suite 309, Cincinnati, Ohio 45227 www.bryanmedical.net


Available Neonate and Pediatric Sizes


I.D. Length	Stock #	I.D. Length	Stock #	I.D. Length	Stock #	I.D. Length	Stock #
2.5 x 26mm	355-2.5-C26	3.5 x 22mm	355-3.5-C22	4.5 x 36mm	355-4.5-C36	5.5 x 40mm	355-5.5-C40
2.5 x 28mm	355-2.5-C28	3.5 x 25mm	355-3.5-C25	4.5 x 39mm	355-4.5-C39	5.5 x 45mm	355-5.5-C45
2.5 x 30mm	350-2.5	3.5 x 28mm	355-3.5-C28	4.5 x 42mm	355-4.5-C42	5.5 x 50mm	355-5.5-C50
2.5 x 32mm	355-2.5	3.5 x 30mm	355-3.5-C30	4.5 x 45mm	355-4.5-C45	5.5 x 55mm	355-5.5
2.5 x 34mm	355-2.5-C34	3.5 x 32mm	355-3.5-C32	4.5 x 48mm	355-4.5	5.5 x 60mm	355-5.5-C60
2.5 x 37mm	355-2.5-C37	3.5 x 34mm	350-3.5	4.5 x 50mm	355-4.5-C50	5.5 x 65mm	355-5.5-C65
2.5 x 40mm	355-2.5-C40	3.5 x 36mm	355-3.5-C36	4.5 x 52mm	355-4.5-C52	5.5 x 70mm	355-5.5-C70
3.0 x 22mm	355-3.0-C22	3.5 x 38mm	355-3.5-C38	4.5 x 56mm	355-4.5-C56	5.5 x 75mm	355-5.5-C75
3.0 x 24mm	355-3.0-C24	3.5 x 40mm	355-3.5	4.5 x 58mm	355-4.5-C58	6.0 x 50mm	355-6.0-C50
3.0 x 26mm	355-3.0-C26	3.5 x 42mm	355-3.5-C42	4.5 x 60mm	355-4.5-C60	6.0 x 54mm	355-6.0-C54
3.0 x 28mm	355-3.0-C28	3.5 x 45mm	355-3.5-C45	4.5 x 65mm	355-4.5-C65	6.0 x 58mm	355-6.0-C58
3.0 x 30mm	355-3.0-C30	3.5 x 50mm	355-3.5-C50	5.0 x 38mm	355-5.0-C38	6.0 x 62mm	355-6.0
3.0 x 32mm	350-3.0	3.5 x 55mm	355-3.5-C55	5.0 x 44mm	355-5.0-C44	6.0 x 66mm	355-6.0-C66
3.0 x 34mm	355-3.0-C34	3.5 x 60mm	355-3.5-C60	5.0 x 50mm	355-5.0	6.0 x 70mm	355-6.0-C70
3.0 x 36mm	355-3.0	4.0 x 30mm	355-4.0-C30	5.0 x 55mm	355-5.0-C55		
3.0 x 38mm	355-3.0-C38	4.0 x 33mm	355-4.0-C33	5.0 x 60mm	355-5.0-C60		
3.0 x 41mm	355-3.0-C41	4.0 x 36mm	350-4.0	5.0 x 65mm	355-5.0-C65		
3.0 x 44mm	355-3.0-C44	4.0 x 38mm	355-4.0-C38	5.0 x 70mm	355-5.0-C70		
3.0 x 47mm	355-3.0-C47	4.0 x 40mm	355-4.0-C40	5.0 x 75mm	355-5.0-C75		
3.0 x 50mm	355-3.0-C50	4.0 x 44mm	355-4.0				
3.0 x 55mm	355-3.0-C55	4.0 x 47mm	355-4.0-C47				
		4.0 x 50mm	355-4.0-C50				
		4.0 x 52mm	355-4.0-C52				
		4.0 x 55mm	355-4.0-C55				
		4.0 x 60mm	355-4.0-C60				
		4.0 x 65mm	355-4.0-C65				

Note: pink denotes standard **neonate**
 tan denotes standard **pediatric**
 ALL other lengths denote custom size.

Measurements			
Size	(1) ID(mm)	(2) OD(mm)	(3) BA
2.5	2.5	3.6	120
3.0	3.0	4.3	120
3.5	3.5	5.0	120
4.0	4.0	5.6	120
4.5	4.5	6.3	120
5.0	5.0	7.0	120
5.5	5.5	7.6	120
6.0	6.0	8.4	120

(1) ID: inside diameter (clear width) of tube;
 (2) OD: outside diameter of tube; (3) BA: bending angle







REF 350-4.0

Size 4,0
OD 5.6mm
ID 4.0
Length 36mm

LOT N1000029769
4035324011995





REF 355-4.0

Size 4,0
OD 5.6mm
ID 4.0
Length 44mm

LOT N1000029797
4035324011994

Color-coded outer packaging:
 TRACOE *mini* neonates: ■ pink
 TRACOE *mini* pediatric: ■ tan

For quick and easy tracking and identification in the ward.